



Recombinant M-phase inducer phosphatase cdc-25.3 (cdc-25.3)

Product Code	CSB-MP335562CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30634
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MCVDVPCENC IVRNDGLRLK CSECAEGSSK LFPRQNRQHS SAISHISNSS PPTRKRSIDG GYTSGTDSAN TSEIVIKKRL TFSKKSHSTS EIETWNAHLQ VDYHLETVTP SCSTVYQKIT SETLIEIMQK LSQIEFMQKY ILIDCRYDYE YNGGHIKGAQ SLFNPETAAD FFFNKDGSKK INRIPIFYCE YSQKRGPTMA NNLREVDKRL NSNIYPRCDY EEIYLLEGGY KNFYAFTRGL EKEQRVQLCE PDNYVIMFDD RYKAELRKHQ FHKKNVSKPM KKWSSTTSVI SILTTSGTRI STLRQTC DPI HEHDAH
Source	Mammalian cell
Target Names	cdc-25.3
Protein Names	Recommended name: M-phase inducer phosphatase cdc-25.3 EC= 3.1.3.48 Alternative name(s): Cell division cycle-related protein 25.3
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.