



Recombinant Putative RNA polymerase II subunit B1 CTD phosphatase R08D7.2 (R08D7.2)

Product Code	CSB-EP335563CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30641
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	METEIPVNEK ETSFRRRVYT TIETLIDIEN TDQLQEHLPH LHCLGWDEVV EERYVNKQCG FPSCQKAPPK ITRNQMFEID RKEGKIFEFR KQRAKFCSEM CYQKSSFVRK QLDEHPLWIT GLTEARTQKV YEVPDES FVS PIAEKSEPVD SQFKKEPSIW LVTDSIIAKV QDMKLSEAE EPKSLDPEDQ DIEPFKLTDD DKDFIKSIKE FRNSNFGPPS TSKLLKTAPK PVLSAKDRKK EDEVLAKLRA KYGKKNALQK KPPILIEAQE IHSKLTMEK AKEAWLVDLI KSWFTPETRK LVREGARPTG GAAEQILMDF LSGKKVDAEK LVNLPNLDKY NVKEKRLNIF LHSIRNHWMD LEARLHLTPT RRDILSRVAS TFQLDSENIT GWTKREINSI VIALFIVICL VDVELGDDYF KKNASPELT AISNELCGLD SFQITGLHAA IKSQC
Source	E.coli
Target Names	R08D7.2
Protein Names	Recommended name: Putative RNA polymerase II subunit B1 CTD phosphatase R08D7.2 EC= 3.1.3.16 Alternative name(s): RNA polymerase II-associated protein 2
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.