



Recombinant Uncharacterized protein ZK643.5 (ZK643.5)

Product Code	CSB-YP335564CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30651
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTTVEEDLLL NEPVATEYDG LTAIDEDFLL GAGDVSRKGT QUESTDDLLGS SEPFQEDIL DTVPEQTDEN EDEAGDDELE SEKEELDYDE EEDDEDRRER TSRYTSEKKG SRKDSVEGDE NKKENGQDET KRSGIPSLFD KTITNNPMKS GDQPIGVILR IEGEAERVFY PPPSFMEPKL NIVPRLPNGI PLIGVPSSHG RAGFGANYQD QARMAAGGGM AGMPMGWENQ VNAFLNNTTK SKSRGRDTRK RRSSSYSSTS SSSDGRSRSR SSSRSDRRRR DDRKRPREDR DRKDYRRDGR RDGRRDDRKR HDDRRRDSYH RDSKDIKRR KEQSKHSAME SAKALGLSND YIDQVNEQKR KREEIVRKKE ERRHAPVSEK KEVPTTVSTN TSSAAYANSK DKTKAYLAVN VTGVQQLPTA VKKIEAIASE LGPIKKCWRS DEDVVSIIFN AHDKAKDFML KHNGKVL SGL RITVSLEKKF LNLNDVN
Source	Yeast
Target Names	ZK643.5
Protein Names	Recommended name: Uncharacterized protein ZK643.5
Expression Region	1-487
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.