



Recombinant Tryptophanase 2 (tnaA2)

Product Code	CSB-BP335593SRW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31015
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	PKGEPFKIK MVEPIRLIPR EDREAAIKAA HYNPFLLRSS DVYIDLLTDS GTGAMSQFQW SAMMLGDESY AGASSYYRLK EAVTDITGYE YVLPTHQGRG AEKSAFAQLI TRPGMYVLSN MFFDTRRGHV QLAGGRPIDL LLDVPTEEYH PFKGNMDTAR LEAFIQEHGA ENIACIVMTV TNNSAGGQPV SMANIRETSR IARKYGILL FDVARYAENC HFIRMREEGY ADKAPIDIAR EMFSYGDGLM MSAKKDALVN IGGLLAFKDE ELYTRVGGTV VPFEGLTYG GLAGRDLEAM AVGLREALDP DYLAIRVGQV EYLGNNLLRSA GIPIQWPVGG HAVFIDAAKF LPHIPWDQFP GHALTVALYQ EGGVRTVEVG SLVMGRDPET GENVRSPFEF TRLAIPRRVY TNLHLEDVAE TVINAFQKRE QIRGVKFTRE PKVLRHFTAH FDLV
Source	Baculovirus
Target Names	tnaA2
Protein Names	Recommended name: Tryptophanase 2 EC= 4.1.99.1 Alternative name(s): L-tryptophan indole-lyase 2 Short name= TNase 2
Expression Region	2-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.