



Recombinant *Saccharomyces cerevisiae* Homoserine dehydrogenase (HOM6)

Product Code	CSB-MP335605SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31116
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSTKVVNVAV IGAGVVGSFAF LDQLLAMKST ITYNLVLLAE AERSLISKDF SPLNVGSDWK AALAASTTKT LPLDDLIAHL KTSPKPVILV DNTSSAYIAG FYTKFVENGI SIATPNKKAF SSDLATWKAL FSNKPTNGFV YHEATVGAGL PIISFLREII QTGDEVEKIE GIFSGTLSYI FNEFSTSQAN DVKFSDDVVKV AKKLGYTEPD PRDDLNGLDV ARKVTIVGRI SGVEVESPTS FVPQSLIPKP LESVKSADF LEKLSDYDKD LTQLKKEAAT ENKVLRFIHK VDVATKSVSV GIEKYDYSHF FASLKGSDNV ISIKTKRYTN PVVIQGAGAG AAVTAAGVLG DVIKIAQRL
Source	Mammalian cell
Target Names	HOM6
Protein Names	Recommended name: Homoserine dehydrogenase Short name= HDH EC= 1.1.1.3
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.