



Recombinant *Saccharomyces cerevisiae* Putative prephenate dehydratase (PHA2)

Product Code	CSB-YP335729SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32452
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MASKTLRVLF LGPKGTYSHQ AALQQFQSTS DVEYLPAASI PQCFNQLEND TSIDYSVVPL ENSTNGQVVF SYDLLRDRMI KKALSLPAPA DTNRITPDIE VIAEQYVPIT HCLISPIQLP NGIASLGNFE EVIIHSHPVQV WGQVECYLRS MAEKFPQVTF IRLDCSSTSE SVNQCIRSST ADCDNILHLA IASETAAQLH KAYIIIEHSIN DKLGNTRFL VLKRRENAGD NEVEDTGLLR VNLFTTTRQ DDPGLVDVL NILKIHSLNM CSINSRPFHL DEHDRNWRYL FFIEYYTEKN TPKNKEKFYE DISDKSKQWC LWGTFPRNER YYHK
Source	Yeast
Target Names	PHA2
Protein Names	Recommended name: Putative prephenate dehydratase Short name= PDT EC= 4.2.1.51 Alternative name(s): Phenylalanine-requiring protein 2
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.