



Recombinant *Saccharomyces cerevisiae* Protein DOM34 (DOM34)

Product Code	CSB-EP335817SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33309
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKVISLKKDS FNKGGAVITL LPEDKEDLFT VYQIVDKDDE LIFKKKFTSK LDEAGKKKST DLVCLKIKVI SEDFDMKDEY LKYKGVTVTD ESGASNV DIP VGKYLSTLD YVYPFTIIKQ NFNKFMQKLL NEACNIEYKS DTAAVVLQEG IAHVCLVTSS STILKQKIEY SMPKKKRTTD VLKFDEKTEK FYKAIYSAMK KDLNFDKLT IILCSPGFYA KILMDKIFQY AEEEHNKKIL DNKGMFFIAH CSTGYLQGIN EVLKNPLYAS KLQDTKYSKE IMVMDEFLLH LNKDDDKAWY GEKEVVKAAE YGAISYLLLT DKVLHSDNIA QREEYKLMD SVESNGGKAL VLSTLHSLGE ELDQLTGIAC ILKYPLPDL EDDGEE
Source	E.coli
Target Names	DOM34
Protein Names	Recommended name: Protein DOM34 EC= 3.1.-.-
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.