



Recombinant Escherichia coli Putative multidrug resistance outer membrane protein MdtQ (mdtQ)

Product Code	CSB-MP335823ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33369
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	CAPMHETRQ ALSQQTPAAQ VDTALPTALK MVGQTXQWWL EYHDNQLTSL INNALQNAPD MQVAEQRIQL AEAQAKAVAT QDGPQIDFSA DMERQKMSAE GLMGPFALND PAAGTTGPWY TNGTFGLTAG WHLDIWGKNR AEVTARLGTV KARAAEREQT RQLLAGSVAR LYWEWQTQAA LNTVLQQIEK EQNTIIATDR QLYQNGITSS VEGVETDINA SKTRQQLNDV AGKMKIIEAR LSALTNNQTK SLKLPVALP KVASQLPDEL GYSLLARRAD LQAAHWYVES SLSTIDAAKA AFYPDINLMA FLQQDALHLS DLFHSAQQM GVTAGLTLPI FDSGRLNANL DIAKAESNLS IASYNKAVVE AVNDVARAAS QVQTLAEKNQ HQAQIERDAL RVVGLAQARF NAGIIAGSRV SEARIPALRE RANGLLLQGQ WLDASIQLTG ALGGGYKR
Source	Mammalian cell
Target Names	mdtQ
Protein Names	Recommended name: Putative multidrug resistance outer membrane protein MdtQ
Expression Region	22-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.