



Recombinant Caprine arthritis encephalitis virus Pol polyprotein (pol)

Product Code	CSB-MP335832CDG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33459
Product Type	Recombinant Protein
Immunogen Species	Caprine arthritis encephalitis virus (strain Cork) (CAEV-Co)
Purity	>85% (SDS-PAGE)
Sequence	TRNHMSQLWK ERTYAKRMQR KERHKGKTAG KREEGDTCGA VRSSYGITSA PPMVQVRIGS QQRNLLFDTG ADRTIVRWHE GSGNPAGRIK LQGIGGIVEG EKWNNVELEY KGETRKGTV VLPQSPVEVL GRDNMARFGI KIIMANLEEK RI
Source	Mammalian cell
Target Names	pol
Protein Names	Recommended name: Pol polyprotein Cleaved into the following 3 chains: 1. Protease Alternative name(s): Retropepsin EC= 3.4.23.- Reverse transcriptase/ribonuclease H Short name= RT EC= 2.7.7.49 EC= 3.1.26.13 Alt
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.