



Recombinant Escherichia coli Uncharacterized protein yejA (yejA)

Product Code	CSB-MP335880ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33913
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	Q AIKESYAFAV LGEPRYAFNF NHFDYVNPAA PKGGQITLSA LGTFDNFNRY ALRGNPGART EQLYDTLFTT SDDEPGSYYP LIAESARYAD DYSWVEVAIN PRARFHDGSP ITARDVEFTF QKFMTEGVPQ FRLVYKGTTV KAIAPLTVRI ELAKPGKEDM LSLFSLPVFP EKYWKDHLKLS DPLATPPLAS GPYRVTSWKM GQNIVYSRVK DYWAANLPVN RGRWNFDTIR YDYYLDDNVA FEFKAGAFD LRMENDAKNW ATRYTGKNFD KKYIHKDEQK NESAQDTRWL AFNIQRPVFS DRRVREAITL ADFEWMNKA LFYNAWSRTN SYFQNTTEYAA RNYPDAAELV LLAPMKKDLP SEVFTQIYQP PVSKGDGYDR DNLLKADKLL NEAGWVLKGQ QRVNATTGQP LSFELLPAS SNSQWVLPFQ HSLQRLGINM DIRKVDNSQI TNRMRSRDYD MPRVWRAMP WPSSDLQISW SSEYINSTYN APGVQSPVID SLINQIIAAQ GNKEKLLPLG RALDRVLTWN YYMLPMWYMA EDRLAWWDFK SQPAVRPIYS LGIDTWWYDV NKAACLPSAS KQGE
Source	Mammalian cell
Target Names	yejA
Protein Names	Recommended name: Uncharacterized protein yejA
Expression Region	20-604
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.