



Recombinant *Saccharomyces cerevisiae* U3 small nucleolar RNA-associated protein 11 (UTP11)

Product Code	CSB-EP335909SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34247
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MAKLVHDVQK KQHRERSQLT SRSRYGFLEK HKDYVKRAQD FHRKQSTLKV LREKAKERNP DEYYHAMHSR KTDAGLLIS SRHGDEEDES LSMDQVKLLK TQDSNYVRTL RQIELKKLEK GAKQLMFKSS GNHTIFVDSR EKMNEFTPEK FFNTTSEMVN RSENRLTKDQ LAQDISNNRN ASSIMPKESL DKKKLLKFKQ VKQHLQRETQ LKQVQQRMDA QRELLKKGSK KKIVDSSGKI SFKWKKQRKR
Source	E.coli
Target Names	UTP11
Protein Names	Recommended name: U3 small nucleolar RNA-associated protein 11 Short name= U3 snoRNA-associated protein 11 Alternative name(s): U three protein 11
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.