



Recombinant Uncharacterized protein B0303.7 (B0303.7)

Product Code	CSB-EP335911CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34258
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MVRLPPPTPP RKPNICGTSG NSNGINENID TSTFSVSAAK ALFESAKRPE VRNAAIAAAK NPFVRQTAKN VAQDEKARGA VLNAVKDPSG GNKISAAFAV ADANRKSDSV PPPPHRGGT SPGKLSGSHS AIRSQLENMH IGG SVMSSNS STATPSYQSN STSSYSSSAP VAPPPVPRHT QSIPSYSAPA ISTPYGIAKF DYAPTQSDM GLRIGDTVLI SKKVDAEWFY GENQNQRFTG IVPSSYLDIK IPLKEAFTAL PPRPAAPSGS SSGTYATAIY DYN SNEAGDL NFAVGSQIMV TARVNEEWLE GECFGRSGIF PSQFVDCPNL YQVPMKQSAP SSAPSYSHP TNNSGPKQTV TVSYDYDSGV ASDLRLFEGD VITALEDDIDA QWLLAE CRGQ QGMVPKTFLG PPYGGSPKKAP SKGIAE IACA FSPKKAVATG DYHSEDPKHL YVTRGDHLLI VEDVDDYYK GKLEAFK TLP AGILPKNIVK LEN
Source	E.coli
Target Names	B0303.7
Protein Names	Recommended name: Uncharacterized protein B0303.7
Expression Region	1-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.