



# Recombinant Uncharacterized protein B0303.7 (B0303.7)

<b>Product Code</b>	CSB-YP335911CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34258
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVRLPPPTTP RKPNICGTSG NSNGINENID TSTFSVSAAK ALFESAKRPE VRNAAIAAAK NPFVRQTAKN VAQDEKARGA VLNAVKDPSG GNKISAAFAV ADANRKSDSV PPPPPHRGGT SPGKLSGSHS AIRSQLENMH IGG SVMSSNS STATPSYQSN STSSYSSSAP VAPPPVPRHT QSIPSYSAPA ISTPYGIAKF DYAPTQSDM GLRIGDTVLI SKKVDAEWFY GENQNQRFTG IVPSSYLDIK IPLKEAFTAL PPRPAAPSGS SSGTYATAIY DYN SNEAGDL NFAVGSQIMV TARVNEEWLE GECFGRSGIF PSQFVDCPNL YQVPMKQSAP SSAPSYSHP TNNSGPKQTV TVSYDYDSGV ASDLRLFEGD VITALEDDA QWLLAECRGQ QGMVPKFTLG PPYGGSPKKAP SKGIAEIACA FSPKKAVATG DYHSEDPKHL YVTRGDHLLI VEDVDDYYK GKLEAFKTLP AGILPKNIVK LEN
<b>Source</b>	Yeast
<b>Target Names</b>	B0303.7
<b>Protein Names</b>	Recommended name: Uncharacterized protein B0303.7
<b>Expression Region</b>	1-493
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.