



Recombinant CID domain-containing protein 1 (cids-1)

Product Code	CSB-EP335915CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34281
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MADFTEQTLR QKLTNLSNHP SSIQTTSAWL LQNHSNRELI IRVWLKTVKK ETHGSKIVNL LYVANDVSQN ARKTCPPQFKD EFFPAIESSF RHAIELKNAK EVEHAIGKLI NVWKDRQIFT PSQCKRLHEV HQQVKLSGSF PTPAVANKEH GKVAQPSQFV VEEAKKNAQD VLLSLKRLQN PPSTEREIRT ELSKYPDNIS CPEKLQSVRS SQEAQSLLVQ NEEALPMLLE YVKRLKNETN ERETLESNLN MLIENVRMSI EHHEKLCREV KRREDRIKAD LLEVEKTFES LPDLAAEMPN APLPTLEALF EKRKK
Source	E.coli
Target Names	cids-1
Protein Names	Recommended name: CID domain-containing protein 1
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.