



Recombinant Uncharacterized protein C05B5.4 (C05B5.4)

Product Code	CSB-MP335917CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34292
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MLRFFQMTWT TALYFLVFLS FGNSMPQKSK PFWLSADINS IEWQEGCLTT VLC SHPRFQL LKDLLPISER ASISWPVTEQ FLEHTVAPFV SYWPSGRIED VLSAQVVGV DTTYGFPRTC DQTPAVRIFP VDLYELAPES AENKTIHLKA KCFEATITVT KHVERCPWCP DPQEILISNE IPQQINLQQS ALESGIHSFV GIFYKSSTSE SLQVGICLLA ALAFFASLAF IIMLSVYLKS KKSSRRNVSV NVQPRILIPYN SRCDDYDTRY DMPWDQQRPL TYWMSKSTVT SPNSLRSEAY QTYRIPPPPN FAPPV
Source	Mammalian cell
Target Names	C05B5.4
Protein Names	Recommended name: Uncharacterized protein C05B5.4
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.