



# Recombinant SET domain-containing protein 3 (set-3)

<b>Product Code</b>	CSB-BP335921CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34318
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVEVPSFRSN EHFQQLCDKV TVKWDKKRGR FVEAIEDIPI GTVVCVETGI TVNVDPQNCY RCLKITENAS FAYCKNCEEF YEPDEIACGE FDELGIFKLA AHLVFSYPFA DIASLVQSSD PEPPSCAPKA LSTQDIEAIF QLTPFPEIGE AFKAPAIQNA IKRIVESLET DENWGRLDQI SRTMTFTKAL RIMAERSAKN AHTIYSIEQI ESQEDNLPMA TGLFPISIF NHSCTPNISG FFVRNTFIFV SQGVRAREEL LDSYGVTYHQ HTFEQRTNFL ASVSGFICHC ESCFKMKSLK VLEKPKKYPE TIATNASCFT DITSYTENIP RGSKELENLI IAFARRPDSE TYTELIFQFW KKFVENAKFR GITYDPYLVR PYVENTILAW NKEVGCENEE KMSLLIVTHR LLRNCYVDLH PLSELVVKLV NVVANQNADE INRTLEKCLK RARMLWKYSE QSEDL
<b>Source</b>	Baculovirus
<b>Target Names</b>	set-3
<b>Protein Names</b>	Recommended name: SET domain-containing protein 3
<b>Expression Region</b>	1-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.