



Recombinant SET domain-containing protein 3 (set-3)

Product Code	CSB-YP335921CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34318
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MVEVPSFRSN EHFPQLCDKV TVKWDKKRGR FVEAIEDIPI GTVVCVETGI TVNVDPQNCY RCLKITENAS FAYCKNCEEF YEPDEIACGE FDELGIFKLA AHLVFSYPFA DIASLVQSSD PEPPSCAPKA LSTQDIEAIF QLTPFPEIGE AFKAPAIQNA IKRIVESLET DENWGRLDQI SRTMTFTKAL RIMAERSAKN AHTIYSIEQI ESQEDNLPMA TGLFPISSIF NHSCTPNISG FFVRNTFIFV SQGVRAREEL LDSYGVTYHQ HTFEQRTNFL ASVSGFICHC ESCFKMKSLK VLEKPKKYPE TIATNASCFT DITSYTENIP RGSKELENLI IAFARRPDSE TYTELIFQFW KKFVENAKFR GITYDPYLVR PYVENTILAW NKEVGCENEE KMSLLIVTHR LLRNCYVDLH PLSELVVKLV NVVANQNADE INRTLEKLKL RARMLWKYSE QSEDL
Source	Yeast
Target Names	set-3
Protein Names	Recommended name: SET domain-containing protein 3
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.