



Recombinant Prion-like- (Q/N-rich) domain-bearing protein 38

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-YP335937CXY |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P34428 |
| Product Type | Recombinant Protein |
| Immunogen Species | Caenorhabditis elegans |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNPNQPPNYP PGNPNQPPPF PPTPYYPMYQ AGMPMPYIPQ SPQYPQPGPP GGPPPFMGIQ IKQEPMWSQE NRYGMQAGPS GMPGMPPSSH PYGYPPMHGM PPGMPQGMPP MMDPNMPDFS SPPMHHTTS LLAQHLSSAP LPPAPPQQQS HLSQLLQQT TYSTATNPHA APTQTIARRP DNEARNAQLQ EQELAAARQQQ QQEYMLQQHQ QQQQNRAGQL QAEQAFANPP TVSIPPQTGH PPQQTSAASSQ QRAQETMAQQ TEPQSHIRES PPKETNLEKY STGELCLYGR ELVNDLNVKT SYLSTILKKV MERKPLNQGE NPNDLADQCK NALQRMSDIR QIIKRREPT WKRMTGEDYI ELMLDDSELK KPMDEERQKR RAELEQRFAI VISKTPPVEG HIWYQGKWIS EELHKKYMQF EANKEIVKNL AAELKALAWK IDVTSPGHLK KVDQTKISKK KTDDQV |
| Source | Yeast |
| Target Names | mdt-30 |
| Protein Names | Recommended name: Prion-like-(Q/N-rich) domain-bearing protein 38 Alternative name(s): Glutamine/asparagine-rich protein pqn-38 |
| Expression Region | 1-466 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.