



Recombinant Putative peroxisomal biogenesis factor 19 (prx-19)

Product Code	CSB-YP335941CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34453
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTDETTQNIK DKTEEELAAL LDQTLGEFTA TPAPKPRTTD DELDELMASA DQEAQAQKAAK DFQKMLEQMV TLQEEAMKKA GADPSEGEGE QPLDPNDPEA LAMMDALKQL MECSSNVANA SNPEEFMAGL DMLRSPNSPM EPFMSMIMQT LASKEVMYPP LKEIFDNYPK YLEDNGAGLD AETKERYEKQ FEVLGKICTE FEKQPELAEV QPVDAAATQPA PEADPASIEH FEKLGKLLVE LQQYGYPPKE LVGALPDGWQ IDESGLPKVA DAAAATEAC
Source	Yeast
Target Names	prx-19
Protein Names	Recommended name: Putative peroxisomal biogenesis factor 19 Alternative name(s): Peroxin-19
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.