



Recombinant Uncharacterized protein F59B2.8 (F59B2.8)

Product Code	CSB-MP335945CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34483
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MEFMSSPFPI TKLPLVPRCK ILKFFDYGDL LDISLCSKRM AQTVRDIHIT ADLHYLTLGK DNQQSIQIQF KQEQKVIYWD FKLDLVEEEM PERRTVGAIV FEKCRKYSQR EIGLFQFSSY HYDIEDSIGE VSKHLFYIFP TSIEIRLSVQ FCRTLRLNLS YEHLQNIDVL TLLGSIMLQN VLEKIFSRVK IRRKLWVEPE TDDEYPILQA LHVENLHLAS ARWMTRDHL QLTCCSVDFH EHYFGWKDLT VFAENWLNNK NSKLEMVRFG WPNDIDLSA NFEGLKTHKW DPKQREKKFL YSKDNRLHRI DCTKGLELER DDGELATWIY EPDPLTSSLY FLVWTERFPE KKRLEKLPKM LAPYYEQLVQ LRKEYDDSCN LEWLLSNPNL KLDEFVETYR ILRGMDSEVR LSSIGRVRRR RIFDEMYNII DAQND
Source	Mammalian cell
Target Names	F59B2.8
Protein Names	Recommended name: Uncharacterized protein F59B2.8
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.