



Recombinant Cell death-related nuclease 6 (crn-6)

Product Code	CSB-MP335949CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34508
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	AYQ CKDNNGSNVD WFVFKLPHL WNHPDNVPIS NGTGFLYFDV NNKNWKLMPQ GMDVENNAVY YTLQQYYNSN MNTTFSYMYN DEWPDSTIWS NSSGHAKGVT VFDQYTGFWM IHSIPKFPSK DMFRFPSNAH YYGQMGICIS YNTVSLATIA QQLFYNTFT YQFNLPQSFA NQFPVLSQLK NKEYNKSPPL TSTKVLKSLG GQHFRHFAKT GEWGKDLYSD FVGPTLKSSI KVETWNHQSG DEYNLPSVCD PNHVQSTMSA KYIRLPY AID YSSYEDHSKF VVAYSESSSK PPIPYVCIGD INRQSHQIHR GGGTMCYDQ ETYFQFANII SETVPCTKAT AEKATLTVLL IAIITFFK
Source	Mammalian cell
Target Names	crn-6
Protein Names	Recommended name: Cell death-related nuclease 6 EC= 3.1.-.-
Expression Region	18-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.