



Recombinant UPF0183 protein T01G9.2 (T01G9.2)

| | |
|--------------------------|---|
| Product Code | CSB-BP335971CXY |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P34692 |
| Product Type | Recombinant Protein |
| Immunogen Species | Caenorhabditis elegans |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNGKQQVTKK PSTEEDLSFT KFKVVPDVGL KSSQFEFVLG MPINQCIAMI QQHPRMLTKV ELKYSKKDPF YQDIIIIYIGS TGIRLYFDGL SQLIKLIEVD NLSMITLTYN DTIFSDPNM ATLDRVNEFF GSTHPGSYDD KHNIYVQSWP GLSFCFPYGG ENSNLEVRPG FGGNLRSLKY DANSQPCLTK MSYRGPNS EPESVDTPFS CYCGQNRTRK VEAIWENGNI VGIDIQFDTQ NGRIVDGEYD VSTYTRQIYF GDSVSDVQSI LGAPTQVFYK SDDKMKIHRG LHKETLYGPP NFFFNYFVMG LDILDFVSK RVVKFVLHTN APGHCDFGMY SRCNFSIFLN DKQYEIRTDK KFDEFSHAFM NDSNPPRPVV LARQEQQPFG STFCYGIKQI IVERTENGFL TSVTIYDGAK EK |
| Source | Baculovirus |
| Target Names | T01G9.2 |
| Protein Names | Recommended name: UPF0183 protein T01G9.2 |
| Expression Region | 1-422 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |