



Recombinant UPF0183 protein T01G9.2 (T01G9.2)

Product Code	CSB-EP335971CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34692
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNGKQQVTKK PSTEEDLSFT KFKVVPDVGL KSSQFEFVLG MPINQCIAMI QQHPRMLTKV ELKYSKKDPF YQDIIIYIGS TGIRLYFDGL SQLIKLIEVD NLSMITLTYN DTIFSDPNM ATLDRVNEFF GSTHPGSYDD KHNIYVQSWP GLSFCFPYGG ENSNLEVRPG FGGNLRSLKY DANSQPCLTK MSYRGPNS EPESVDTPFS CYCGQNRTRK VEAIWENGNI VGIDIQFDTQ NGRIVDGEYD VSTYTRQIYF GDSVSDVQSI LGAPTQVFYK SDDKMKIHRG LHKETLYGPP NFFFNYFVMG LDILDFVSK RVVKFVLHTN APGHCDFGMY SRCNFSIFLN DKQYEIRTDK KFDEFSHAFM NDSNPPRPVV LARQEQQPFG STFCYGIKQI IVERTENGFL TSVTIYDGAK EK
Source	E.coli
Target Names	T01G9.2
Protein Names	Recommended name: UPF0183 protein T01G9.2
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.