



Recombinant *Saccharomyces cerevisiae* Protein KES1 (KES1)

Product Code	CSB-BP336093SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35844
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSQYASSSSW TSFLKSIASF NGDLSSLSAP PFILSPISLT EFSQYWAEHP ELFLEPSFIN DDNYKEHCLI DPEVESPELA RMLAVTKWFI STLKSQYCSR NESLGSEKKP LNPFLGELFV GKWENKEHPE FGETVLLSEQ VSHHPPVTAF SIFNDKNKVK LQGYNQIKAS FTKSLMLTVK QFGHTMLDIK DESYLVTPPP LHIEGILVAS PFVELEGKSY IQSSTGLLCV IEFSGRGYFS GKKNSFKARI YKDSKDSKDK EKALYTISGQ WSGSSKIKA NKKEESRLFY DAARIPAEHL NVKPLEEQHP LESRKAWYDV AGAIKLGDFN LIAKTKTELE ETQRELKREE EAKGISWQRR WFKDFDYSVT PEEGALVPEK DDTFLKLASA LNLSTKNAPS GTLVGDKEDR KEDLSSIHWR FQRELWDEEK EIVL
Source	Baculovirus
Target Names	KES1
Protein Names	Recommended name: Protein KES1 Alternative name(s): Oxysterol-binding protein homolog 4
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.