



Recombinant Enterobacteria phage P1 Recombination enhancement function protein (ref)

Product Code	CSB-BP336105EDD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35926
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage P1 (Bacteriophage P1)
Purity	≥85% (SDS-PAGE)
Sequence	MKTIEQKIEQ CRKWQKAARE RAIARQREKL ADPVWRESQY QKMRDTLDRR IAKQKERPPA SKTRKSAVKI KSRGLKGRTP TAEERRIANA LGALPCIACY MHGVISNEVS LHIIAGRTAP GCHKKQLPLC RWHHQHAAPA EVREKYPWLW PVHADGVVGG KKEFTLLNKS EMELLADAYE MANIMH
Source	Baculovirus
Target Names	ref
Protein Names	Recommended name: Recombination enhancement function protein
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.