



Recombinant *Saccharomyces cerevisiae* DOA4-independent degradation protein 4 (DID4)

Product Code	CSB-BP336128SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36108
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSLFEWVFGK NVTPQERLKK NQRALERTQR ELEREKRLKLE LQDKKLVSEI KKSANNGQVA AAKVQAKDLV RTRNYIQKFD NMKAQLQAIS LRIQAVRSSD QMTRSMSEAT GLLAGMNRTM NLPQLQRISM EFEKQSDLMG QRQEFMDEAI DNVMGDEVDE DEEADEIVNK VLDEIGVDLN SQLQSTPQNL VSNAPIAETA MGIPEPIGAG SEFHGNPDDD LQARLNTLKK QT
Source	Baculovirus
Target Names	DID4
Protein Names	Recommended name: DOA4-independent degradation protein 4 Alternative name(s): ESCRT-III complex subunit VPS2 Vacuolar protein-sorting-associated protein 2 Vacuolar protein-targeting protein 14
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.