



Recombinant *Saccharomyces cerevisiae* Protein PET10 (PET10)

Product Code	CSB-BP336135SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36139
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SESSISSK PSVELPQATW SHLQRYPALS KFIKYAESLP PVERLISFNL VVLGSVNQWV SESSSSPRLV KQVVAAGKEG AFKLDELVNL LVFKEGVDGL LYNWKSHSNT PGIWLVWFFV DYVANISNTL LREFLIKPLH LQGSTASKEI GSSGEENKVT DASSLPHVAE LSSTTRGMSQ EIQSKVKSNI IDPTKDLAKE KYDAIVKPTT DKLQSVYIDP TKTKLNETYQ RFTTVYENNL SKSESVPKAI VSTGLDLGNA TIEKLGASRE DQTNSKPAAV STN
Source	Baculovirus
Target Names	PET10
Protein Names	Recommended name: Protein PET10 Alternative name(s): Petite colonies protein 10
Expression Region	2-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.