



Recombinant *Saccharomyces cerevisiae* DASH complex subunit DAD2 (DAD2)

Product Code	CSB-BP336139SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36162
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDSIDEQIAI KRKELQSLQK ITSLTDGLKI QLTELNEQIK EMGMNADSV QLMNNWDSII NNISQASLGL LQYAEGDYEI GPWKDSKKKE SEQSNETGLE AQENDKNDDED NDEDEDLVPL PETMVRIRVD GNE
Source	Baculovirus
Target Names	DAD2
Protein Names	Recommended name: DASH complex subunit DAD2 Alternative name(s): DUO1 and DAM1-interacting protein 2 Helper of ASK1 protein 1 Outer kinetochore protein DAD2
Expression Region	1-133
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.