



Recombinant Schizosaccharomyces pombe 26S protease regulatory subunit 4 homolog (mts2)

Product Code	CSB-MP336193SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36612
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGQAQSGNFS NFGDGANGDN KKDQKKDKPK YEPPVPTRTG RRKKKAQSGP DASAKLPTVI PTTRCRLRL KMQRIHDHLL MEEYVQNQE RLKPQDERTQ EERNRVDEIR GTPMSVGTLE EIIDDDHAIV STAGPEYYVS IMSFVDKML EPGCSVLLHH KAMSIVGLLL DDTDPMINVM KLDKAPTESY ADIGGLESQI QEIKEAVELP LTHPELYEEM GIKPPKGVIL YGAPGTGKTL LAKAVANQTS ATFLRVVGSE LIQKYLGDGP RLVRQLFNAA EEHSPSIVFI DEIDAIGTKR YDAQSGAERE IQRTMLELLN QLDGFDTSQR DIKVMATNR ISDLDPALIR PGRIDRKILF ENPDEATKRK IFTIHTSKMN LGEDVNLEEL IQCKDDLGA EIKAIVSEAG LLALRERRMR VVMDDFRQAR EKVLKTKDEG GPAGGLYI
Source	Mammalian cell
Target Names	mts2
Protein Names	Recommended name: 26S protease regulatory subunit 4 homolog Alternative name(s): Protein mts2
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.