



Recombinant Escherichia coli Ribonucleoside-diphosphate reductase 2 subunit beta (nrdF)

Product Code	CSB-BP336265ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37146
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKLSRISAIN WNKISDDKDL EVWNRLTSNF WLPEKVPLSN DIPAWQTLTV VEQQLTMRVF TGLTLLDTLQ NVIGAPSLMP DALTPHEEAV LSNISFMEAV HARSYSSIFS TLCQTKDVDA AYAWSEENAP LQRKAQIIQQ HYRGDDPLKK KIASVFLESF LFYSGFWLPM YFSSRGKLTN TADLIRLIIR DEAVHGYIIG YKYQKNMEKI SLGQREELKS FAFDLLLELY DNELQYTDEL YAETPWADDV KAFLCYNANK ALMNLGYEPL FPAEMAEVNP AILAALSPNA DENHDFFSGS GSSYVMGKAV ETEDEDWNF
Source	Baculovirus
Target Names	nrdF
Protein Names	Recommended name: Ribonucleoside-diphosphate reductase 2 subunit beta EC= 1.17.4.1 Alternative name(s): R2F protein Ribonucleotide reductase 2
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.