



Recombinant Solanum tuberosum Cytochrome b-c1 complex subunit Rieske, mitochondrial (FES1)

Product Code	CSB-EP336358FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37841
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	≥85% (SDS-PAGE)
Sequence	SSNSVSP AHQLGLVSDL PATVAAIKNP SSKIVYDDSN HERYPPGDPS KRAFAYFVLT GGRFVYASSV RLLILKFVLS MSASKDVLAL ASLEVDLSSI EPGSTVTVKW RGKPVFIRRR TDDIKLANS VDLGTLRDPQ QDAERVKNPE WLVVVGVCTH LGCIPLNAG DFGGWFCPCH GSHYDISGRI RKGPPAPYNLE VPTYSFLEEN KLLIG
Source	E.coli
Target Names	FES1
Protein Names	Recommended name: Cytochrome b-c1 complex subunit Rieske, mitochondrial EC= 1.10.2.2 Alternative name(s): Complex III subunit 5 Rieske iron-sulfur protein Short name= RISP Ubiquinol-cytochrome c reductase iron-sulfur subunit
Expression Region	54-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.