



Recombinant Bacillus subtilis HTH-type transcriptional regulator degA (degA)

Product Code	CSB-EP336373BRJ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37947
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	>85% (SDS-PAGE)
Sequence	MKTTIYDVAK AAGVSITTVS RVINNTGRIS DKTRQKVMNV MNEMAYTPNV HAAALTGKRT NMIALVAPDI SNPFYGELAK SIEERADELG FQMLICSTDY DPKKETKYFS VLKQKKVDGI IFATGIESHD SMSALEEIAS EQIPIAMISQ DKPLLPMDIV VIDDVRGGYE AAKHLLSLGH TNIACIIGDG STTGEKNRIK GFRQAMEEAG VPIDESLIQ TRFSLESGKE EAGKLLDRNA PTAIFAFNDV LACAAIQAAR IRGIKVPDDL SIIGFDNTIL AEMAAPPLTT VAQPIKEMGA ERHRTAGRSN RGKRKAKQKI VLPPELVVRH STSPLNT
Source	E.coli
Target Names	degA
Protein Names	Recommended name: HTH-type transcriptional regulator degA Alternative name(s): Degradation activator
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.