



Recombinant *Saccharomyces cerevisiae* E3 SUMO-protein ligase MMS21 (MMS21)

Product Code	CSB-BP336471SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38632
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MALNDNPIPK SVPLHPKSGK YFHNLHARDL SNIYQQCYKQ IDETINQLVD STSPSTIGIE EQVADITSTY KLLSTYESES NSFDEHIKDL KKNFKQSSDA CPQIDLSTWD KYRTGELTAP KLSELYLNMP TPEPATMVNN TDTLKILKVL PYIWNDPTCV IPDLQNPAD EDDLQIEGGKI ELTCPITCKP YEAPLISRKC NHVFDRDGIQ NYLQGYTTRD CPQAACSQVV SMRDFVRDPI MELRCKIAKM KESQEQDKRS SQAIDVL
Source	Baculovirus
Target Names	MMS21
Protein Names	Recommended name: E3 SUMO-protein ligase MMS21 EC= 6.3.2.- Alternative name(s): Methyl methanesulfonate-sensitivity protein 21 Non-structural maintenance of chromosome element 2 Short name= Non-SMC element 2
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.