



Recombinant *Saccharomyces cerevisiae* Proteasome-interacting protein CIC1 (CIC1)

Product Code	CSB-EP336492SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38779
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAKKSNSKKS TPVSTPSKEK KKVIEKKSST AIPRERVIKA VNELIKFTSK PQDENNEEGN NGKKNLLEDD EEELKKDLQL IVVNNKSFTG TSKSFKLKL NVKHSFYKPW KEASATAVKD FKVLLILKDS DIKKVSEDDL FDQLDSEGIK VDEIICGKDL KTVYKAYEAR NAFISQFSLI LADDSIVTSL PKLMGGKAYN KVETTPISIR THANKEFSLT TLTNNIKKVY MNQLPVKLPR GTTLNVHLGN LEWLRPEEFV DNVELISEQL IKAYQIRSIF IKTNRSPVLP LYYNQDVLDE LEAKKDKIEE THEDDMVTID GVQVHLSTFN KGLMEIANPS ELGSIFSKQI NNAKKRSSSE LEKESSESEA VKKAKS
Source	E.coli
Target Names	CIC1
Protein Names	Recommended name: Proteasome-interacting protein CIC1 Alternative name(s): Core interacting component 1
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.