



Recombinant *Saccharomyces cerevisiae* Protein LSM12 (LSM12)

Product Code	CSB-EP336499SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38828
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	SVSLEQTLG FRIKVTNVLD VVTEGRLYSF NSSNNTLTIQ TTKKNQSPQN FKVIKCTFIK HLEVIGDKPS FNSFKKQQIK PSYVNVERVE KLLKESVIAS KKKELLRGKG VSAEGQFIFD QIFKTIGDTK WVAKDIILD DVKVQPPYKV EDIKVLHEGS NQSITLIQRI VERSWEQLEQ DDGRKGG
Source	E.coli
Target Names	LSM12
Protein Names	Recommended name: Protein LSM12 Alternative name(s): Sm-like protein 12
Expression Region	2-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.