



Recombinant *Saccharomyces cerevisiae* Protein DSE2 (DSE2)

Product Code	CSB-YP336501SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38844
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	NSNGETVKL ITSDGIVYSY AVYTKTLAPA RVVVKTISYT TTRVYPITLA NSVVSSTTEK ITEVSTVSAS EQVSATQTNS LVSTSTVSTI SPTISSGSST SSSSTYDIES SQSIESSGTS SATAEPSASS GFRLTSSSAF VSSTAPFSSQ LSSSSSSETS SSSFSTSSSS APLSLTSSSS SSSSFATIIT LAPSSSKSGN SQLTLASSSS TSAVESSQTG STIARTTSTL VPSSSVDTTS RATTSMPLS SSTQSISSVSS SDGTCYVFYD DDDYYSTVYL TNPSQSVD
Source	Yeast
Target Names	DSE2
Protein Names	Recommended name: Protein DSE2 Alternative name(s): Daughter-specific expression protein 2
Expression Region	22-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.