



# Recombinant *Saccharomyces cerevisiae* FK506-binding nuclear protein (FPR3)

<b>Product Code</b>	CSB-BP336511SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P38911
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SDLLPLATY SLNVEPYTPV PAIDVTMPIT VRITMAALNP EAIDEENKPS TLRIIKRNPD FEDDDFLGGD FDEDEIDEES SEEEEEKTQ KKKKSKGKKA ESESEDDEED DDEDDEFQES VLLTLSPEAQ YQQLDLTIT PEEEVQFIVT GSYAISLSGN YVKHPFDTPM GVEGEDEDED ADIYDSEDYD LTPDEDEIIG DDMDDLDEE EEEVRIEEVQ EEDEEDNDGE EEQEEEEEEE QKEEVKPEPK KSKKEKKRKH EEKKEEKKAK KVKKVEFKKD LEEGPTKPKS KKEQDKHKPK SKVLEGGIVI EDRTIGDGPQ AKRGARVGM R YIGKLKNGKV FDKNTSGKPF AFKLG RGEVI KGWDIGVAGM SVGGERRIII POPYAYGKQA LPGIPANSEL TFDVKLVSMK N
<b>Source</b>	Baculovirus
<b>Target Names</b>	FPR3
<b>Protein Names</b>	Recommended name: FK506-binding nuclear protein EC= 5.2.1.8 Alternative name(s): FKBP-70 Nucleolar proline isomerase Peptidyl-prolyl cis-trans isomerase Short name= PPIase Proline rotamase
<b>Expression Region</b>	2-411
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.