




# Recombinant Arabidopsis thaliana Nucleoside diphosphate kinase 1 (NDK1)

<b>Product Code</b>	CSB-YP336542DOA
<b>Relevance</b>	Major role in the synthesis of nucleoside triphosphates other than ATP. The ATP gamma phosphate is transferred to the NDP beta phosphate via a ping-pong mechanism, using a phosphorylated active-site intermediate. Plays a role in response to reactive oxygen species (ROS) stress.
<b>Abbreviation</b>	Recombinant Mouse-ear cress NDK1 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P39207
<b>Alias</b>	Nucleoside diphosphate kinase I Short name:NDK I Short name:NDP kinase I Short name:NDPK I
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MEQTFIMIKPDGVQRGLIGEVICRFEKKGFTLKGLKLISVERSFAEKHYEDLSSK SFFSGLVDYIVSGPVVAMIWEGKNVVLTRKIIIGATNPAASEPGTIRGDFIDIG RNVIHGSDSVESARKEIALWFPDGPVNWQSSVHPWVYET
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	NDK1
<b>Protein Names</b>	Nucleoside diphosphate kinase I Short name:NDK I Short name:NDP kinase I Short name:NDPK I
<b>Expression Region</b>	1-149aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	20.5 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.