



Recombinant Escherichia coli Uncharacterized KamA family protein YjeK (yjeK)

Product Code	CSB-YP336553ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39280
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MAHIVTLNTP SREDWLTQLA DVVTDPELL RLLNIDAE EK LLAGRS AKKL FALRVPR SFI DRMEKGNPDD PLLRQVLTSQ DEFVIAPGFS TDPLEEQHSV VPGLLHKYHN RALLLVKGGC AVNCRYCFRR HFPYAENQGN KRNWQTALEY VAAHPELDEM IFSGGDPLMA KDHELDWLLT QLEAIPHIKR LRIHSRLPIV IPARITEALV ECFARSTLQI LLVNHINHAN EVDETFRQAM AKLRRVGVTL LNQSVLLRDV NDNAQTLANL SNALFDAGVM PYYLHVLDKV QGAAHFMVSD DEARQIMREL LTLVSGYLVP KLAREIGGEP SKTPLDLQLR QQ
Source	Yeast
Target Names	epmB
Protein Names	Recommended name: Uncharacterized KamA family protein YjeK
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.