



Recombinant Escherichia coli Uncharacterized protein yjgR (yjgR)

Product Code	CSB-MP336559ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39342
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSEPLLIART PDTELFLLPG MANRHGLITG ATGTGKTVTL QKLAESLSEI GVPVFMADV K GDLTGVAQAG TVSEKLLARL KNIGVNDWQP HANPVVVWDI FGEKGHVPVRA TVSDLGPLLL ARLLNLNDVQ SGVLNIIIFRI ADDQGLLLLL FKDLRAITQY IGDNAKSFQN QYGNISSASV GAIQRGLLSL EQQGAAHFFG EPMLDIKDWM RTDANGKGI NILSAEKLYQ MPKLYAASLL WMLSELYEQL PEAGDLEKPK LVFFFDEAHL LFNDAPQVLL DKIEQVIRLI RSKGVGVWVWV SQNPSDIPDN VLGQLGNRVQ HALRAFTPKD QKAVKAAAQT MRANPAFDTE KAIQELGTGE ALISFLDAKG SPSVVERAMV IAPCSRMPGV TEDERNGLIN HSPVYGKYED EVDRESAYEM LQKGFQASTE QQNNPPAKGK EVAVDDGILG GLKDILFGTT GPRGGKKGDV VQTMAKSAAR QVTNQIVRGM LGSLLGGRRR
Source	Mammalian cell
Target Names	yjgR
Protein Names	Recommended name: Uncharacterized protein yjgR
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.