



Recombinant Escherichia coli Probable HTH-type transcriptional regulator YjjM (yjjM)

Product Code	CSB-MP336569ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39399
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MSRSQNLNRHN VINQVIDDMA RGHIPSPLPS QSALAEMYNI SRTTVRHILS HLRECGVLTQ VGNDYVIARK PDHDDGFACT TASMSEQNKV FEQAFFTMIN QRQLRPGETF SELQLARAAG VSPVVVREYL LKFGRYNLIH SEKRGQWSMK QFDQSYAEQL FELREMLETH SLQHFLNLPD HDPRWLQAKT MLERHLLRD NIGNSFRMFS QLDRDFHSL LSAADNIFFD QSLEISVIF HFHYQWDESD LKQRNIIAVD EHMTILSALI CRSDLDTLA LRNHLNSAKQ SMIRSINENT RYAH
Source	Mammalian cell
Target Names	IgoR
Protein Names	Recommended name: Probable HTH-type transcriptional regulator YjjM
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.