



Recombinant Enterobacteria phage T4 Uncharacterized 12.1 kDa protein in dda-modA intergenic region (y01A)

Product Code	CSB-EP336572EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39419
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MVYVYAIVYR DKDGFTAPVP LDEHRPAVFF EWKIADKVFT TLKEQYRLAL GKGIPRLVET PRKFWFNKIE VKHVKPDVDT QRLYQRILDT GRIVSIPIAG NLR
Source	E.coli
Target Names	y01A
Protein Names	Recommended name: Uncharacterized 12.1 kDa protein in dda-modA intergenic region
Expression Region	1-103
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.