



Recombinant Guinea pig Alkyldihydroxyacetonephosphate synthase, peroxisomal (AGPS), partial

| | |
|--------------------------|---|
| Product Code | CSB-YP001457GU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P97275 |
| Product Type | Recombinant Protein |
| Immunogen Species | Cavia porcellus (Guinea pig) |
| Purity | >85% (SDS-PAGE) |
| Source | Yeast |
| Target Names | AGPS |
| Protein Names | Recommended name: Alkyldihydroxyacetonephosphate synthase, peroxisomal Short name= Alkyl-DHAP synthase EC= 2.5.1.26 Alternative name(s): Alkylglycerone-phosphate synthase |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Target Details | This gene is a member of the FAD-binding oxidoreductase/transferase type 4 family. It encodes a protein that catalyzes the second step of ether lipid biosynthesis in which acyl-dihydroxyacetonephosphate (DHAP) is converted to alkyl-DHAP by the addition of a long chain alcohol and the removal of a long-chain acid anion. The protein is localized to the inner aspect of the peroxisomal membrane and requires FAD as a cofactor. Mutations in this gene have been associated with rhizomelic chondrodysplasia punctata, type 3 and Zellweger syndrome. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |