



# Recombinant *Saccharomyces cerevisiae* PHO85 cyclin-6 (PCL6)

<b>Product Code</b>	CSB-YP336666SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P40038
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSIKGDSPSS TNASSSPKST YSIQSDDKAN LGSGNVDIRT DNSQQDSNNR RDIVVVTRVA SEETLESQSS TSSMGIRPES SFNYEDASNQ ARVEMNNRVH GSNMNTINKY YPVRFPKNE RQLSDTNNLN EKVQGTHTVQ SSTQEDKILD GDTSNSQVTP SLNIAEFPTD KLLKMLTALL TKIISNDRT AATNPSTLQE IENGRCLALS DNEKKYLSPV LGFRGKHVPQ IGLDQYFQRI QKYCPTTNDV FLSLLVYFDR ISKRCNSVTT TPKTNTAKHE SPSNESSLDK ANRGADKMSA CNSNENNEND DSDDENTGVQ RDSRAHPQMF VMDSHNIHRL IIAGITVSTK FLSDFYNSNS RYSRVGGISL QELNHLELQF LVLCDFELLI SVNELQRYAD LLYRFWNNAK AQSQUALVTGM
<b>Source</b>	Yeast
<b>Target Names</b>	PCL6
<b>Protein Names</b>	Recommended name: PHO85 cyclin-6
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.