



# Recombinant Aspergillus fumigatus Heat shock protein 90 (hsp90), partial

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP336701NGS   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P40292  |
| <b>Storage Buffer</b>    | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0  |
| <b>Product Type</b>      | Recombinant Proteins  |
| <b>Immunogen Species</b> | Aspergillus fumigatus (strain ATCC MYA-4609 / CBS 101355 / FGSC A1100 / Af293) (Neosartorya fumigata)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | hsp90   |
| <b>Protein Names</b>     | Recommended name: Heat shock protein 90<br>Alternative name(s): 65 kDa IgE-binding protein Heat shock protein hsp1<br>Allergen= Asp f 12  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Partial   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |