



Recombinant Drosophila melanogaster S-adenosylmethionine synthase (Sam-S)

Product Code	CSB-BP336704DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40320
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPQKTNGHSA NGCNGSNGNS YDMEDGATFL FTSESVGEGH PDKMCDQISD AILDAHLKQD PNAKVACETV AKTGMILLCG EITSKAVVDY QKVVRETVQH IGYDDSSKGF DFKTCNVLLA LDQQSPEIAA GVHVNRAEEE IGAGDQGIMF GYATDETEEC MPLTVVLAHK LNEKIAELRR SDVFWWARPD SKTQVTCEYL FNQGSVAVPKR VHTIVVSMQH SEKISLET LR SEVMEKVVKV VIPAKYIDAN TIVHINPCGL FVIGGPMGDA GLTGRKIIVD TYGGWGAHGG GAFSGKDFTK VDRSAAYAAR WVAKSLVKAG LCKRCLVQVS YAIGLAEPLS ITVFDYGTSH KSQKELLDII RRNFDLRPGK IVKDLNLRQP IYQRTSTYGH FGRAGFSWEE AKPLEIDN
Source	Baculovirus
Target Names	Sam-S
Protein Names	Recommended name: S-adenosylmethionine synthase Short name= AdoMet synthase EC= 2.5.1.6 Alternative name(s): Methionine adenosyltransferase Short name= MAT
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.