



Recombinant *Schizosaccharomyces pombe* Meiotic recombination protein rec12 (rec12)

Product Code	CSB-BP336714SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40384
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNSNDKKKVV RSWIEQFVHD FVEQLSKPTK DSVNVALKRR KHNSWNGSLD SKANERQKVK VFSFPRNETT IAQLFRVLDV VHEAVISDTV ITKRDIYYRD VDLFRQRTVV DELLGDISNT IGCSRSDLNV EASAKGLVFG SIHIALENGT VITATKPLLI SHHRISITS TAKWVLVIEK EAVFQTLTEE ALADTIIVTA KGFPLMTRK FLVKLAKALP DAKFFGIFDW DPHGLCIYSC FKYGSNAYSH EPHSQLRNLQ LLGPLYEDIF NKNQEFSLKL NKRDIKMITT LLQFEGFQKE PVVREQLQRM LFIQKKAIEIQ AILEFPSWIK GKLADADKSG KHSVR
Source	Baculovirus
Target Names	rec12
Protein Names	Recommended name: Meiotic recombination protein rec12
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.