



Recombinant *Drosophila melanogaster* Myosin regulatory light chain sqh (sqh)

Product Code	CSB-MP336720DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40423
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSSRKTAGRR ATTKKRAQRA TSNVFAMFDQ AQIAEFKEAF NMIDQNRDGF VEKEDLHDML ASLGKNPTDD YLDGMMNEAP GPINFTMFLT LFGERLQGTD PEDVIKNAFG CFDEENMGVL PEDRLRELLT TMGDRFTDED VDEMYREAPI KNGLFDYLEF TRILKHGAKD KDEQ
Source	Mammalian cell
Target Names	sqh
Protein Names	Recommended name: Myosin regulatory light chain sqh Alternative name(s): Myosin regulatory light chain, nonmuscle Short name= MRLC-C Protein spaghetti-squash
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.