



Recombinant *Saccharomyces cerevisiae* Protein SIM1 (SIM1)

Product Code	CSB-MP336728SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40472
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	L PHVDVHQEDA HQHKRAVAYK YVYETVVVDS DGHTVTPAAS EVATAATSAI ITTSVLAPTS SAAAADSSAS IAVSSAALAK NEKISDAAAS ATASTSQGAS SSSSSSSATS TLESSSVSSS SEEAAPTSTV VSTSSATQSS ASSATKSSTS STSPSTSTST STSSTSSSSS SSSSSSSSSS GSGSIYGDLA DFSGPSEKFQ DGTIPCDKFP SGQGVISIDW IGEWWGSGVE NTDSTGGSC KEGSYCSYSC QPGMSKTQWP SDQPSDGRSV GLLCKNGYL YRSNTDADYL CEWGVAAAYV VSKLSKGVAI CRTDYPGTEN MVIPTYVEGG SSLPLTVVDQ DTYFTWEGKK TSAQYYVNNA GSVEDGCIW GTSGSGIGNW APLNFGAGST GGVTYLSLIP NPNNSDALNY NVKIVAADDS SNVIGECVYE NGEFSGGADG CTVSVTSGKA HFVLYN
Source	Mammalian cell
Target Names	SIM1
Protein Names	Recommended name: Protein SIM1
Expression Region	20-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.