



# Recombinant Salmonella typhimurium Secreted effector protein sopD (sopD)

<b>Product Code</b>	CSB-EP336770SXB-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P40722
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPVTLSFGNH QNYTLNESRL AHLLSADKEK AIHMGGWDKV QDHFRAEKKD HALEVLHSII HGQGRGEPGE MEVNVEDINK IYAFKRLQHL ACPAHQDLFT IKMDASQTQF LLMVGDTVIS QSNIKDILNI SDDAVIESMS REERQLFLQI CEVIGSKMTW HPELLQESIS TLRKEVTGNA QIKTAVYEMM RPAEAPDHPL VEWQDSLAD EKSMLACINA GNFEPTTQFC KIGYQEVQGE VAFSMMHPCI SYLLHSYSPF SEFKPTNSGF LKKLNQDYND YHAKKMFIDV ILEKLYLTHE RSLHIGKDGC SRNILLT
<b>Source</b>	E.coli
<b>Target Names</b>	sopD
<b>Protein Names</b>	Recommended name: Secreted effector protein sopD Alternative name(s): Salmonella outer protein D
<b>Expression Region</b>	1-317
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.